

Transferrin Mouse Monoclonal Antibody(1B12)

Catalog No: RA10064

Basic Information	
Host species	Mouse
Applications	IHC,ELISA
Species Cross-Reactivity	Н
Specificity	The Transferrin Mouse Monoclonal Antibody can detects endogenous Human Transferrin protein.
Recommended dilutions	IHC: 1:100-200 ELISA:1:100,00-50,000
	Optimal dilutions should be determined by the end user.
Applications	
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal
Background	
Alternative Names	Apotransferrin antibody, beta 1 metal binding globulin antibody, serotransferrin
	antibody, TF antibody
Observed band	77
Human Gene ID	7018
Human Swiss-Prot Number	P02787
	Transferrins are iron-binding blood plasma glycoproteins that control the level of
Background	free iron in biological fluids. It is responsible for the transport of iron from sites of absorption and heme degradation to those of storage and utilization.