

PPAR Delta Mouse Monoclonal Antibody(2F9)

Catalog No: RA10080

Basic Information

Host species	Mouse
Applications	WB,IHC
Species Cross-Reactivity	H,R,M
Specificity	The PPAR Delta antibody can detects endogenous PPAR Delta protein.
Recommended dilutions	WB: 1:1,000-2,000 IHC: 1:100-200 Optimal dilutions should be determined by the end user.


Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal


Background

Alternative Names	FAAR, NR1C2, NUC1, Peroxisome proliferative activated receptor delta
Observed band	50
Human Gene ID	5467
Human Swiss-Prot Number	Q03181
Background	.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using PPAR Delta Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using PPAR Delta Mouse mAb diluted at 1:200.