

# PPAR Delta Mouse Monoclonal Antibody(1D7)

## Catalog No: RA10082

### Basic Information

Host species	Mouse
Applications	IHC
Species Cross-Reactivity	H,R,M
Specificity	The PPAR Delta antibody can detects endogenous PPAR Delta protein.
Recommended dilutions	IHC: 1:100-200 Optimal dilutions should be determined by the end user.


### Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal


### Background

Alternative Names	FAAR, NR1C2, NUC1, Peroxisome proliferative activated receptor delta
Observed band	50
Human Gene ID	5467
Human Swiss-Prot Number	Q03181
Background	.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using PPAR Delta Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using PPAR Delta Mouse mAb diluted at 1:200.