


## **Caspase-3 Mouse Monoclonal Antibody(4B8)**

Catalog No: RA10105


Basic Information	
Host species	Mouse
Applications	IHC
Species Cross-Reactivity	H,R,M
Specificity	Antibody can detects endogenous Active Caspase 3 protein.
Recommended dilutions	IHC: 1:100-200
	Optimal dilutions should be determined by the end user.
Applications	
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal
Background	
Alternative Names	APOPAIN antibody, CASP 3 antibody, CPP32 antibody, Cysterine protease CPP32 antibody
Observed band	17、32
Human Gene ID	836
Human Swiss-Prot Number	P42574
	Caspase-3 (CPP-32, Apoptain, Yama, SCA-1) is a critical executioner of apoptosis,
Background	as it is either partially or totally responsible for the proteolytic cleavage of many
	key proteins such as the nuclear enzyme poly(ADP-ribose) polymerase (PARP).



## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using Caspase-3 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse LiverTissue using Caspase-3 Mouse mAb diluted at 1:200.