

Muscle Actin Mouse Monoclonal Antibody(2H3)

Catalog No: RA10127

Basic Information

Host species	Mouse
Applications	IHC
Species Cross-Reactivity	H,R,M
Specificity	The Muscle Actin antibody can detects endogenous Muscle Actin protein.
Recommended dilutions	IHC: 1:100-200 Optimal dilutions should be determined by the end user.


Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal


Background

Alternative Names	
Observed band	42
Human Gene ID	58
Human Swiss-Prot Number	P68133
Background	

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded Mouse Skeletal Muscle Tissue using Muscle Actin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Muscle Actin Mouse mAb diluted at 1:200.