


## HP-1 $\gamma$ Rabbit Polyclonal Antibody (HP-1 g)

Catalog No: RA20280


Basic Information	
Host species	Rabbit
Applications	WB,IHC
Species Cross-Reactivity	H,R,M
Specificity	Antibody can detects endogenous HP-1γproteins.
Recommended dilutions	WB: 1:200-500 IHC: 1:200-500
	Optimal dilutions should be determined by the end user.
Applications	
Formulation	Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal
Background	
Alternative Names	CBX3, Chromobox homolog 3, HECH, HP1 gamma
Observed band	24
Human Gene ID	11335
Human Swiss-Prot Number	Q13185
Background	Heterochromatin protein 1 (HP1) is a family of heterochromatic adaptor molecules involved in both gene silencing and higher order chromatin structure. All three HP1 family members ( $\alpha$ , $\beta$ , and $\gamma$ ) are primarily associated with centromeric heterochromatin.



## Selected Validation Data



Western blot analysis of 1) Hela Cell Lysate, 2) Rat Brain Tissue Lysate using HP-1y Rabbit pAb diluted at 1:1000.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using HP-1 $\gamma$  Rabbit pAb diluted at 1:500.