


ATG7 Rabbit Polyclonal Antibody (F309)

Catalog No: RA20297


Basic Information	
Host species	Rabbit
Applications	IHC
Species Cross-Reactivity	H, R, M
Specificity	The Antibody can detects endogenous ATG7proteins.
Recommended dilutions	IHC: 1:100-200
	Optimal dilutions should be determined by the end user.
Applications	
Formulation	Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% sodium azide as
	Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal
Background	
Alternative Names	APG7, GSA7, Autophagy related protein 7
Observed band	78
Human Gene ID	10533
Human Swiss-Prot Number	O95352
Background	The molecular machinery of autophagy was largely discovered in yeast and referred to as autophagy-related (Atg) genes. Formation of the autophagosome involves a ubiquitin-like conjugation system in which Atg12 is covalently bound to Atg5 and targeted to autophagosome vesicles. This conjugation reaction is mediated by the ubiquitin E1-like enzyme Atg7 and the E2-like enzyme Atg10.



Selected Validation Data



Immunohistochemical analysis of paraffin-embedded Human Hepatocarcinoma Tissue using ATG7 Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using ATG7 Rabbit pAb diluted at 1:200.