

# GAL4 Activation Domain Rabbit Polyclonal Antibody(AD)

**Catalog No: RA20378**

## Basic Information

Host species	Rabbit
Applications	WB
Species Cross-Reactivity	N/A
Specificity	Antibody can detects recombinant GAL4 Activation Domain protein.
Recommended dilutions	WB: 1:10,000-50,000 Optimal dilutions should be determined by the end user.

## Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal

## Background

Alternative Names	N/A
Observed band	N/A
Human Gene ID	
Human Swiss-Prot Number	N/A
Background	GAL4 is a transcriptional activator which controls the expression of genes involved in the catabolism of galactose.

Selected Validation Data



Western blot analysis of Recombinant AD protein with GAL4  
Activation Domain Rabbit pAb diluted at 1:20,000.