

## **ECFP Rabbit Polyclonal Antibody**

Catalog No: RA20405

Basic Information	
Host species	Rabbit
Applications	WB, IP
Species Cross-Reactivity	N/A
Specificity	The ECFP Antibody can detects ECFP & CFP tag protein
Recommended dilutions	WB: 1:2,000-5,000 IP: 1:200
	Optimal dilutions should be determined by the end user.
Applications	
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal
Background	
Alternative Names	Cyan fluorescent protein, Enhanced Cyan fluorescent protein
Observed band	N/A
Human Gene ID	
Human Swiss-Prot Number	N/A
Background	Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish Aequorea victoria.



## Selected Validation Data



Western blot analysis of ECFP recombinant protein with ECFP Rabbit Polyclonal Antibody diluted at 1:5,000.