

Goat anti-Rabbit IgG (H+L),HRP Conjugated

Catalog No: RSA10

Basic Information

Host species	Goat
Applications	WB,IHC,ELISA
Species Cross-Reactivity	Rabbit
Recommended dilutions	WB: 1:4,000-12,000 IHC: 1:500-2,500 ELISA:1:4,000-20,000 Optimal dilutions should be determined by the end user.

Applications

Formulation	0.01M PBS, 50% glycerol, pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	2 mg/ml
Clonality	Polyclonal
Note	Please do not use in dirty and bright environments to avoid pollution and failure. In addition, this reagent is only used for scientific research and cannot be directly used for diagnosis.