

## Goat anti-Mouse IgG (H+L),FITC Conjugated

Catalog No: RSA13

### Basic Information

Host species	Goat
Applications	IHC, IF/FC
Species Cross-Reactivity	Mouse
Recommended dilutions	IHC: 1:500-2,500 FC/ IF:1:150-1,500 Optimal dilutions should be determined by the end user.

### Applications

Formulation	0.01M PBS, 0.5% BSA, pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal
Note	Please do not use in dirty and bright environments to avoid pollution and failure. In addition, this reagent is only used for scientific research and cannot be directly used for diagnosis.