

His-Tag Rabbit Polyclonal Antibody


Catalog No: RTA54

Basic Information	
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	N/A
Specificity	The His tag antibody can recognize C-terminal, internal, and N-terminal His-tag fusion proteins.
Recommended dilutions	WB: 1:5,000
	Optimal dilutions should be determined by the end user.
Applications	
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal
Background	
Alternative Names	
Observed band	N/A
Human Gene ID	N/A
Human Swiss-Prot Number	N/A
Background	A polyhistidine-tag is an amino acid motif in proteins that consists of at least five histidine (His) residues, often at the N- or C-terminus of the protein. Monoclonal antibodies specific to six histidine tags can greatly improve the effectiveness of several different kinds of immunoassays, helping researchers identify, detect, and purify polyhistidine fusion proteins in bacteria, insect cells, and mammalian cells.



Antibodies · ELISA kits · Biochemical reagents · Recombinant proteins · Immune related reagents

Selected Validation Data



1
ug His fusion protein+ Primary antibody dilution at 1
 $_{\rm 1:5,000}$ 2 $_{\rm 1:10,000}$