

RFP-Tag Rabbit Polyclonal Antibody

Catalog No: RTA58

Basic Information

Host species	Rabbit
Applications	WB
Species Cross-Reactivity	N/A
Specificity	The RFP tag antibody can recognize RFP tag fusion proteins.
Recommended dilutions	WB: 1:5,000 Optimal dilutions should be determined by the end user.

Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal

Background

Alternative Names	
Observed band	N/A
Human Gene ID	N/A
Human Swiss-Prot Number	N/A
Background	The red fluorescent protein cloned from Discosoma coral (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the Aequorea jellyfish.

Selected Validation Data

