

## SBP-tag Rabbit Polyclonal Antibody

## Catalog No: RTA60

Basic Information	
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Recombinant
Specificity	Antibody can detect endogenous Recombinant SBP protein.
Recommended dilutions	WB: 1:5,000-10,000
	Optimal dilutions should be determined by the end user.
Applications	
Formulation	Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal
Background	
Alternative Names	Streptavidin Binding Peptide, CEQLRARLEHHPQGQREP
Observed band	ALL
Human Gene ID	
Human Swiss-Prot Number	
Background	Streptavidin Binding Peptide (SBP) tag is very versatile and can be inserted into expression vectors for easy detection of recombinant proteins and purification by use of affinity columns.



Antibodies · ELISA kits · Biochemical reagents · Recombinant proteins · Immune related reagents

## Selected Validation Data



Western blot analysis of Recombinant SBP-Tag protein using SBP -Tag Rabbit pAb diluted at 1:20,000.

www.reedbiotech.com